Absolut-Drucksensor U11320

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10/08 Hh

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10/08 Hh

1. Safety instructions

- To avoid permanent damage to the built-in semiconductor sensor, never exceed the maximum threshold pressure of 1000 kPa.
- Suitable only for non-corrosive gases such as air, helium and nitrogen.
- Do not allow the sensor element to come into contact with water.

2. Description

The absolute pressure sensor with its extensive measurement range is particularly suitable for experiments to demonstrate Boyle's law and for measuring the piston pressure (PV diagram) in a Stirling engine. In addition, the absolute pressure sensor can also be used to record and measure the production of oxygen during photosynthesis and for transpiration experiments in a closed system.

Two-port measurement procedure for the pressure sensor: nozzle 1 is connected to the external pressure via a connecting nipple, and nozzle 2 is connected to a sealed reference vacuum.

The sensor box is automatically recognised via the interface.

3. Equipment supplied

- 1 Sensor box
- 1 8-pin mini DIN connection lead, length: 60 cm
- 1 Silicone tube, inner dia.: 2 mm, length: 1 m
- 1 Plastic syringe 20ml

4. Technical data

Measurement range: 0 to 250 kPa
Sensor type: Semiconductor sensor
Accuracy: $\pm 1\%$
Resolution: 0.1 Pa
Connections: Serrated nozzle 4.8 mm dia.

5. Operation

- Use the full length of the silicone tube or shorten it to the length desired.
- Use the silicone tube to connect the pressure source to the nozzle of the sensor.
- During the experiment, the elasticity of the tube should be taken into account – this could possibly lead to an error in readings.

6. Sample experiment

6.1 Measuring the absolute pressure in relation to the volume (Boyle's law)

Apparatus required:

1	3B NETlog™ interface	U11300
1	Absolute pressure sensor	U11320

- Set-up the experiment according to Fig. 1.
- Fit the plastic syringe with an approx. 2-cm long hose.
- Fill the syringe with 20 ml of ordinary air.

- Push the free end of the silicone tube onto the nozzle of the pressure sensor.
- CAUTION: as far as possible, do not alter the volume in the syringe!
- Connect the absolute pressure sensor to the 3B NETlog™ interface and wait for the interface to recognise the sensor.
- The first reading appears on the interface display.
- In the 3B NETlab™ software's manual mode, enter the readings for the volume by hand in steps of 1 ml at a time.
- Plot the graph of the characteristic.

Fig. 1: Measuring the absolute pressure in relation to the volume


Fig. 2: Pressure against volume

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