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
3/03 ALF

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
Instruction sheet

3/03 ALF

- ① Spring clip
- ② Base track
- ③ Glass rod
- ④ Test rods
- ⑤ Pointer
- ⑥ Scale

The apparatus demonstrates the linear expansion of solid materials and allows the determination of expansion coefficients for copper, iron and glass.

1. Safety instructions

- Caution. The experimental procedure results in hot steam.
- Do not touch heated rods with your fingers. Use cloths when replacing rods.
- Do not subject the glass rod to mechanical stresses.

2. Description, technical data

The equipment consists of a base track with a spring clip at the left-hand end for attaching a testing rod. At a distance of 50 cm from the end there is a notch in the base strip for the pointer. The copper and iron test rods have a ring nut 65 cm from the end for placing on the pointer. The glass rod has a metal ring with a ring nut in the same place. Behind the pointer is a 0-5 cm scale. To feed in hot steam, a 10 cm long glass rod attached to a hose is provided.

| | |
|--------------------|-------------------------|
| Dimensions: | 530 mm x 60 mm x 240 mm |
| Weight: | 0.6 kg |
| Length of rods: | 630 mm approx. |
| Diameter of rod: | 8 mm approx. |
| Length of pointer: | 200 mm |
| Scale markings: | mm |
| Scale factor: | 1 : 50 |

3. Principle

To determine the linear expansion coefficient α for various materials, it is necessary to measure the expansion for a certain temperature rise ΔT . Thus the rods are heated to 100° C by means of steam and the temperature difference ΔT from room temperature is calculated. The expansion is determined from the movement of the pointer d , where a pointer movement of 50 mm indicates an expansion of 1mm. The expansion coefficient is determined from the extension w (scale factor 1:50) and the length of rod l between the two fixed points by means of the formula:

$$\alpha = \frac{d}{l \cdot w \cdot \Delta T}$$

4. Instructions for use

Also required for heating the rods is a vapor generator or a bunsen burner heating an Erlenmeyer flask

- Attach the end of the rod without the nut to a rubber hose and secure in the spring clip.
- Place the pointer in the notch under the scale and attach the upper end of the pointer to the rod with the ring nut.
- Adjust the pointer to zero by sliding the rod.
- Attach to a steam generator or an Erlenmeyer flask half-filled with water by means of the short glass rod and a hose.
- Boil water. Steam flows through the test rod and heats it to approximately 100° C.

(**Note:** at high altitudes, water boils at slightly less than 100° C.)

- When steam has been passing through the rod for about 1 minute and no more condensing steam is emerging from the end of the rod, read off the largest movement of the pointer.

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