

## U20600 Kundt's tube

## U20601 Microphone probe

## U20602 Battery box

### Operating instructions

11/02 ALF



The equipment set comprising Kundt's tube and accessories is meant to display stationary sound waves with open or closed tube ends and determine wavelengths in air and other gases.

#### 1. Safety instructions

- Protect the microphone and loudspeaker against moisture.
- The external voltage through the microphone's connection line should not exceed 5V.
- Do not clean the acrylic glass body with aggressive agents or solvents.

#### 2. Description, technical data

The equipment set designated Kundt's tube consists of an acrylic-glass tube with a scale, two removable end plates and an integrated hose nipple for filling the tube with various gases. One end plate is furnished with a loudspeaker, the other with a bore and guide for


mounting a movable piston or microphone probe (U20601).

The equipment set includes two clamps for mounting Kundt's tube on a tripod, and cables for connecting the loudspeaker.

Length:	1000 mm
Diameter:	70 mm
Hose nipple:	7 mm Ø
Scale:	1000 mm
Division:	mm and cm

#### Drawing:

- ① End plate with loudspeaker, 4-mm jacks and hose nipple
- ② Resonance tube
- ③ End plate with bore and guide for mounting a piston or microphone probe
- ④ Microphone probe
- ⑤ Microphone
- ⑥ Scale



## 2.2. Microphone probe

The microphone probe is used to measure changes in sound pressure inside Kundt's tube.

A miniature microphone is attached to the end of a long rod made of stainless steel. It is connected by means of a 5-pole DIN plug to the battery box (U20602). This battery box also has a terminal for connecting an oscilloscope or voltmeter. The microphone probe can be connected directly to the digital counter (U21000) via the adapter U20603.




Microphone's frequency range: 20 Hz to 20000 Hz  
 Probe's dimensions: 740 mm x 8 mm Ø  
 Connection cable's length: 2 m

## 2.3. Battery box

The battery box supplies the microphones (for instance, U20601 or U18030) and other analog sensors with a power of 5 V DC so that they can be connected directly with a measuring device or an oscilloscope.

The box consists of a compartment for a 9-V alkaline battery which supplies the required 5 V DC via a regulator. Two 6-pole DIN jacks (180°) and one 8-pole DIN jack (270°) are available as input channels. Two 5-pole DIN jacks and two 4-mm safety jacks serve for connecting measuring devices.

Dimensions: 143 mm x 84 mm x 37 mm



- ① DIN jacks for connecting measuring devices
- ② 4-mm safety jacks for connecting an oscilloscope, voltmeter or interface
- ③ DIN jack for connecting various sensors or a LabPro interface via an output adapter
- ④ DIN jacks for connecting microphones


**Note:** In order to make simultaneous use of the microphone probe U20601 and an oscilloscope, connect the microphone probe to the sensor input (3) and the oscilloscope to the output (2).

## 3. Sample experiments


### 3.1 Stationary waves in a closed tube



### 3.2 Stationary waves in carbon dioxide



### 3.3 Stationary waves in a tube with one closed end




### 3. Stationary waves in an open tube



### 3.5 Changes in the air column

Sound source: Tuning fork or loudspeaker



To perform these experiments, additional use is required of a function generator (for example, U21015) to excite the loudspeaker, and an oscilloscope (for instance, U11175) to display the oscillation nodes and antinodes.

### 3.6 Determination of the speed of sound in air

- Install Kundt's tube with closed ends and the microphone probe on a tripod; connect a function generator (for instance, U21015) and an oscilloscope (for instance, U11175).
- Apply a frequency  $f = 2700$  Hz to the loudspeaker.
- Move the microphone probe past the nodes and measure the intervals between them using the scale.
- The average interval between two neighbouring nodes turns out to be 6.3 cm.
- This results in a wavelength  $\lambda = 12.6$  cm.
- The speed of sound  $c$  is calculated with the formula  $c = f \cdot \lambda$   
 $c = 2,7 \cdot 10^3 \cdot 12,6 \cdot 10^{-2} \text{ m/s} = 340 \text{ m/s}$