


## Leslie's Cube U8442835

### Instruction sheet

01/11 ALF



#### 1. Description

The Leslie's cube is used for investigating heat radiation from a hot body as a function of its temperature, colour and surface characteristics.

Leslie's cube is a hollow aluminium cube with a removable lid for pouring in hot water and 2 openings for inserting a thermometer or thermal sensor and a stirrer. One side each is polished, matt, white and black.

#### 2. Scope of delivery

- 1 Leslie's Cube
- 1 Stirrer
- 2 Rubber stoppers with 6 mm aperture

#### 3. Technical data

Dimensions: 100x100x100 mm<sup>3</sup> approx.  
 Weight: 360 g approx.

#### 4. Operation

To conduct the experiment, the following articles are additionally recommended:

- |                            |          |
|----------------------------|----------|
| 1 Thermopile               | U8441301 |
| 1 Multimeter ESCOLA 10     | U8531161 |
| 1 Thermometer              | U14295   |
| or                         |          |
| 1 Digital thermometer      | U11853   |
| and                        |          |
| 1 Immersion sensor NiCr-Ni | U11854   |
| 1 Stand base               | U8611210 |
| 2 Experiment leads         |          |

- Take off the lid and fill the cube with water or oil (not exceeding approx. 130° C).
- Insert a thermometer for monitoring the temperature and the stirrer into the appropriate apertures. Close the lid again.
- The liquid used to warm the cube from within must be well stirred to achieve an even heat distribution.
- Set up the thermopile approx. 3 to 5 cm away from Leslie's cube.
- Connect the multimeter to the thermopile and set the mode switch to the smallest DC current range. (The current measuring range will produce a larger needle deflection in the comparison to the voltage measuring range).
- After setting up the experiment, wait for a few minutes before taking readings.

Note:

Readings may be made incorrect due to the influence of body heat or other external influences.

- Do not touch the thermopile while taking readings.
- Avoid direct sunlight and do not set up the apparatus in the vicinity of a heater/radiator.

The black and white lacquered faces emit approximately equal amounts of infrared radiation. The reason for this is that white and black only appear as such within the visible light spectrum. If only the emitted thermal radiation is observed, which has longer wavelengths than visible light, both the white and black faces appear as so-called grey bodies; in other words, both radiate all wavelengths within this range with equal intensity. By contrast, the thermal emission of the metallic faces is much weaker.




Fig 1 Experimental set-up