


## 8497180 Transformer core with yoke and clamps

## 8497200 Pair of drilled pole pieces

### Instruction Sheet

05/05 JH



- ① Transformer core (U-core)
- ② Yoke (I-core)
- ③ Clamps for clamping the yoke (the bar that completes the ring) or pole pieces (8497200) firmly to the core.
- ④ Pair of drilled pole pieces

The transformer core with yoke and clamps are provided for use in conjunction with the accessories listed under item 4 of the instructions for assembling the demountable transformer.

The pole pieces are required for experiments on electromagnetism where a well-defined air gap is necessary (e.g. Waltenhofen's pendulum or investigations of paramagnetic and diamagnetic samples).

#### 1. Safety instructions

- Polished surfaces should be kept free of dirt and grease.
- Transformer core, yoke and pole pieces should not be exposed to moisture.
- For transporting the U-core and yoke (I-core), make sure the clamps are firmly secured.
- During the experiment, the yoke or pole pieces should be firmly secured using the clamps.

#### 2. Description, technical data

##### 2.1 Transformer core 8497180

Transformer core and yoke made of high-quality, laminated iron for use in transformers, with two holes drilled for securing pole pieces or yoke with the aid of clamps.

Cross section of core: 40 mm x 40 mm  
 Height including yoke: 170 mm  
 Width: 150 mm  
 Material: Laminated iron  
 Weight: 6 kg approx.

##### 2.1.1 Scope of delivery:

Transformer core  
 Yoke  
 Clamps (pair)

## 2.2 Pair of drilled pole pieces 8497200

Pole pieces made of soft iron each with one plane and one conical end. The pole pieces have a bore drilled through the middle of them.

Cross section of core: 40 mm x 40 mm  
 Bore at conical end: 5 mm  
 Bore at plane end: 12 mm  
 Material: Soft iron  
 Weight: 1.7 kg approx.

## 3. Operation

- The safety instructions for the coils must be observed at all times.
- Mount the primary and secondary coils on the core.
- Lay the polished side of the yoke or pole pieces on top of the U-core.
- Attach the clamps.
- Firmly secure the yoke or pole pieces with the clamps.

## 4. Example experiments



Step-up transformer



Waltenhofen's pendulum



Spark discharge along horn-shaped electrodes

### Accessories for step-up transformer

Item	Cat. no.	Winding turns	Tap(s)
Mains coil 220 V	8497420	600	
Mains coil 115 V	Available on request	300	
High current coil for needle-point melting experiments	8497406	6	
Coil	8497410	72	6/30/54/66/72
Coil	8497430	600	200
Coil	8497440	1200	400
Coil	8497450	6000	2000
Melting ladle	8497310	1	
Coil with 5 winding turns High current coil for spot welding	8497320	5	
High-tension coil with horn-shaped spark electrodes	8497460	24000	

## 5. Cleaning, maintenance and storage

- Polished surfaces should be kept free of dirt and grease.
- Store in dry conditions.
- Remove any rust with fine steel wool or sandpaper.