*...going one step further*

**1013460**

## I. Introduction

Most people associate flavour enhancers with the term glutamate. In fact, glutamate or its derivatives are the most important flavour enhancers worldwide. Typically, glutamate is added to pre-cooked meals, frozen foods, snacks, tinned food, salad dressings and meat or fish dishes. In some countries it is even common to put it on the table as a condiment for the seasoning of food.

The Japanese chemist Kikunae Ikeda discovered in 1907 that there was a fifth flavor besides sweet, sour, bitter and salty, which he called umami (from Japanese „umai“ meaty, savoury, tasty). He found that the umami taste has to do with the L-glutamic acid and its salts. He recognized its importance for enhancing the flavour of food, so that two years later, he founded a company for the production and selling of flavour enhancers together with a partner from the industry sector.

In both food and the human body, glutamic acid can occur in two mirror-image isomers, as L-glutamic acid and as D-glutamic acid, wherein the L-form predominates in volume. Furthermore, L-glutamic acid exists in bound and free form. In bound form, glutamic acid is combined with other amino acids and forms proteins. In free form, glutamic acid is present as a single amino acid. Only the free glutamic acid or its salts play an important role for the umami taste of foods. Most foods contain glutamic acid primarily in bound form. However, in soybeans, tomatoes, milk, spinach and many cheeses, for instance, glutamic acid occurs naturally in high concentrations.

The greatest importance with regard to flavour enhancers is attached to monosodium glutamate. Other salts of glutamic acid are monopotassium glutamate, calcium diglutamate, and monoammonium glutamate and magnesium diglutamate.

The detection of these flavour enhancers forms the basis of this experimental kit.

For completeness, it should be mentioned that there are flavour enhancers which are not based on glutamate, but share the same flavour characteristics (umami taste). These are guanylic acid, disodium guanylate, dipotassium guanylate, calcium guanylate, inosinic acid, disodium inosinate, dipotassium inosinate, calcium inosinate, calcium 5'-ribonucleotides and disodium 5'-ribonucleotides.

In Europe, flavour enhancers are additives which need to be declared either by name or with their corresponding E number in the list of ingredients in foods (see table 1). Australia and New Zealand have adopted the European labeling standards. In the U.S. and Canada one can also find similar rules for flavour enhancers.

# Tracking Flavour Enhancers

English

**Table 1:** The E number list which is widely used in Europe. The various flavour enhancers are listed together with their corresponding E number.

E Number	Designation
E620	Glutamic acid
E621	Monosodium glutamate
E622	Monopotassium glutamate
E623	Calcium diglutamate
E624	Monoammonium glutamate
E625	Magnesium diglutamate
E626	Guanylic acid
E627	Disodium guanylate
E628	Dipotassium guanylate
E629	Calcium guanylate
E630	Inosinic acid
E631	Disodium inosinate
E632	Dipotassium inosinate
E633	Calcium inosinate
E634	Calcium 5'-ribonucleotide
E635	Disodium 5'-ribonucleotide
E636	Maltol
E637	Ethylmaltol
E640	Glycine and its sodium salt

E636, E637 and E640 do not cause umami taste.

E636 (maltol) enhances sweet taste.

E636 (ethylmaltol) is similar to maltol, but more intense in its effect.

E640 (glycine) tastes sweet.

Flavour enhancers are controversial. In many studies on humans and animals, negative effects are described. There are also studies which show no causal relationship, however these were mostly launched by the large manufacturers of flavour enhancers.

Nevertheless, with flavour enhancers one can earn much money. Food manufacturers' interest in them is in fact huge because certain expensive ingredients such as meat, shrimp or cheese can be replaced or reduced by using them. Several renowned scientists are very worried as they see correlations between flavour enhancers and the onset and development of Alzheimer's disease, ADHD (attention deficit-hyperactivity disorder), obesity and other diseases.

It should also be mentioned that the addition of glutamate to baby food is totally forbidden in Germany, Sweden and other European countries. In the USA, producers have agreed on a voluntary ban of glutamate in baby food.

Since many consumers now respond more negatively to the addition of flavour enhancers to foods, many

# Tracking Flavour Enhancers

manufacturers stopped using the pure additive (E620 ff) in their products, but use instead a mixture of various amino acids or their salts. In many countries it is sufficient to list these ingredients as “spice”, “yeast extract”, “fermented wheat” or the like on the label. Ultimately all these terms are „Trojan horses“ used to hide the fact that flavour enhancers were added to the product.

## II. Components of the Kit

Extraction solution, 250 ml  
L-glutamic acid, 20 mg  
Ninhydrin, 300 mg  
Filter papers, 20 pieces  
Glass capillaries (volume 10 µl), 50 pieces  
Chromatography foil, 20 pieces  
Chromatography chamber, 5 pieces

## III. Additional Equipment and Solutions Required

Erlenmeyer flasks, 100 ml  
Mortar  
Funnel  
Tweezers  
Bowl (glass or plastic)  
Hair dryer, alternatively oven (105°C)  
H<sub>2</sub>O dist., sterile  
Disposable gloves (powder-free), suitable for laboratory work

## IV. Preparation of Reagents

### L-Glutamic Acid

Dissolve 20 mg glutamic acid in 20 ml sterile H<sub>2</sub>O (dist.). This gives a concentration of 1 mg/ml. This concentration can be further diluted to 1:3 or 1:5 to get standards with 0.33 or 0.2 mg/ml. The glutamic acid standards should be stored frozen at –18°C in aliquots.

### Plasticizer Medium

The plasticizer medium for thin layer chromatography is composed of butanol, acetic acid and H<sub>2</sub>O in the ratio of 4:1:1.

### Ninhydrin Solvent

The ninhydrin solvent (100 ml) is composed of 95 ml 2-propanol and 5 ml acetic acid.

### Preparation of Ninhydrin

Dissolve 300 mg ninhydrin in 100 ml ninhydrin solvent. Please use a clean Erlenmeyer flask. The solution can be stored in a closed bottle in the refrigerator.

## V. Experiment

### Sample Preparation

Weigh 3 g of a sample and put it into an Erlenmeyer flask. Mix it with 5 ml of extraction solution and finely ground it in a mortar. If there is not a large enough mortar available, the sample may also be processed portion by portion. Filter the extract through a paper filter. The filtrate shall be used for the thin-layer chromatography.

If you want to carry out the thin-layer chromatography later, you can store the filtrate frozen at  $-18\text{ }^{\circ}\text{C}$  for several weeks.

### Thin-Layer Chromatography

**Note:** Please do not touch the chromatography paper with your fingers. Please wear powder-free gloves use tweezers and touch only the edge of the paper when moving it. The ruler that you need for later experimental steps should be absolutely clean!

First, draw a thin horizontal line with a pencil about 1.5 cm from the bottom of the chromatography paper (see Fig 1). This is the starting line for the samples and the glutamine acid standard. Below this starting line, at spaces of 8 mm, mark the numbers 1,2,3 etc. from left to right where the traces will run. Part of the above mentioned glutamic acid standard (1 mg/ml) shall be diluted 1:3 and 1:5 diluted with  $\text{H}_2\text{O}$ , so that then 3 standards with the concentrations 1mg/ml, 0.33 mg/ml and 0.2 mg/ml are available. Fill the chromatography chamber with the solvent mixture to about 7 mm and seal it with the lid. For the saturation of the chamber, it should be incubated for at least 15 min.

Using capillary tubes, apply both the samples and the standards (10  $\mu\text{l}$  each) on the starting line at equally spaced intervals (above the numbers). For each sample and each standard a new capillary tube must be used.

For application, place the chromatography paper horizontally on the table, then open the chromatography chamber, and place the chromatography paper loaded with samples and standards with the starting line pointing downwards into the chamber. Close it again with the lid.

Now the solvent mixture rises upward in the chromatography paper. When the solvent has risen about 7cm (measured from the starting line), remove the paper from the chamber and put it on the table to dry. The drying process can be accelerated using a hair dryer on a low setting or in an oven. The chromatogram now ready to be stained, however if you prefer it is possible to store the chromatogram in the refrigerator protected from light for a few days.

### Staining of Amino Acids

**Note:** The bowl has to be absolutely clean (no fingerprints etc.) as any protein contamination would be stained, too!

Pour some ninhydrin staining solution into a suitable bowl (glass or plastic). Then pull the chromatography paper very briefly (1 sec.) through the staining solution and let the paper dry


in the air or an oven at about  $105^{\circ}\text{C}$ . Now the chromatogram is ready for analysis.

The stained chromatogram can be stored at room temperature in the dark. Sticking it in a workbook may best be done with scotch tape.

## VI. Evaluation

On the chromatogram you can see several spots in various shades from red to violet. The glutamic acid standards are shown in lanes 6, 7 and 8 of fig. 1. The glutamate in the samples (here in this aqueous environment as glutamic acid) runs on the same level as the standards, which makes it easily identifiable. By comparing the spot size of the sample with the spot sizes of the standards, the amount of glutamate in the samples can be semi-quantitatively estimated. The spots above and below the third row represent different amino acids that cannot be identified with this kit system.

**Figure 1:** Example of a chromatogram with samples (lanes 1-5) and three glutamic acid standards (lanes 6-8). Concentrations of glutamic acid standard: 0.2 mg/ml, 0.33 mg/ml and 1 mg/ml. The sample in lane 1 contains no glutamic acid, the samples in lanes 2-5 show varying amounts of glutamic acid. The estimation of the amount of glutamic acid shall be carried out in comparison with the standards.

## VII. Suggestions and Questions for the Integration into Teaching

### Preparatory work before the lesson:

- 1 Do some research on the topic „The flavour enhancer glutamate.“ (Alternatively, a lecture on the subject might be prepared and held by a student or teacher.)
- 2 Do some research on the topic of „TLC“. (Alternatively, a lecture on the subject might be prepared and held by a student or teacher.)
- 3 For the next lesson, bring a small amount of food into the classroom, for which the presence of the flavour enhancer glutamate shall be investigated.

### Tasks during the lesson:

- 1 Read the entire experiment. Clarify questions with your neighbour.
- 2 Explain the function of the glutamic acid standard.
- 3 Justify why three glutamic-acid standards of different concentrations are used.
- 4 Describe the experiment's underlying issue.
- 5 Carry out the experiment.
- 6 Log the experiment. Attach the chromatogram with tape to your protocol and label it.
- 7 Evaluate the experiment with regard to the question in task 2.
- 8 Check results you did not expect at all by searching for possible errors in the experimental procedures. Please read the manual again. If necessary, run the experiment once more.

# Tracking Flavour Enhancers

English

## Solutions:

### for task no. 2:

- The glutamic acid serves as a standard control approach. Without it, a glutamic acid spot of the sample mixture can not be identified because there are other amino acids present in the tested food, these would also be stained.

### no. 3:

- Using the dilutions, the amount of glutamic acid in the studied samples can be estimated semi-quantitatively.

### no. 4:

- For which kind of food can you prove the flavour enhancer MSG?
- In which kind of food can you detect the greatest concentration of glutamate, in comparison to other foods examined?

## VIII. Safety and Disposal

The safe handling of chemicals and laboratory equipment requires a certain level of basic knowledge and safety measures. During the experiment, lab coats and protective goggles should be worn. Gloves must be provided by the teacher and should be used if necessary.

Also the handling of the equipment and risks of its use should be known. In particular, one should pay attention to the dangers of electricity. All connectors, power cords, work areas and hands must be dry before handling electrical devices.

Other health and safety measures: tie up long hair, remove any jewellery, and wear tight-fitting sleeves, to avoid unwanted contact with equipment, chemicals, etc.

All waste should be disposed of in accordance with the instructions and local regulations.

## Hazard Classification of the Kit Components

### Extraction Solution

Classification according to Regulation (EC) No 1272/2008:

Not harmful, contains H<sub>2</sub>O and NaCl merely in the physiological range.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC:

Not harmful, contains H<sub>2</sub>O and NaCl merely in the physiological range.

### L-Glutamic Acid

Classification according to Regulation (EC) No 1272/2008:

Not harmful.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC:

Not harmful.

## **Ninhydrin**

Classification according to Regulation (EC) No 1272/2008:

### Hazards

H302: Harmful if swallowed.

H315: Causes skin irritation.

H319: Causes eye irritation.

H335: May cause irritation to the respiratory system.

### Safety instructions

P261: Avoid inhalation

P305 + P351 + P338: On contact with eyes, rinse the eyes gently with water for several minutes. Remove contact lenses if necessary, continue rinsing.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC:

### Hazards

R22: Harmful if swallowed.

R36/37/38: Irritating to eyes, respiratory system and skin.

### Safety instructions

S26: Avoid contact with eyes, rinse immediately with water and seek medical advice.

S28: After contact with skin, wash immediately with plenty of water.

S36: Wear suitable protective clothing.

## **IX. Literature**

Hahn-Deinstorp, E. 2000. Applied Thin-Layer Chromatography – Best Practice and Avoidance of Mistakes. New York: Wiley-VCH.

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© Copyright 2012 for instruction manual and design of product:
**3B Scientific GmbH**
3B Scientific GmbH, Germany