## Relativ-Drucksensor $\pm 1000$ hPa U11323

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03/08 Hh

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03/08 Hh

#### 1. Safety instructions

- To avoid permanent damage to the internal semiconductor sensor, the maximum permitted relative pressure of 4000 hPa must never be exceeded.

Only suitable for use with non-corrosive gases such as air, helium or nitrogen.

- Do not allow the sensor element to come into contact with water.

#### 2. Description

Relative pressure sensor with a measurement range up to 1000 hPa, suitable for measuring the pressure on the piston of the transparent Stirling engine U10050 (for a pV diagram).

For two-port measurement using the sensor, hose connections are provided for two inputs.

The sensor box is designed to be detected automatically by the 3B NETlog™ unit.

#### 3. Equipment supplied

- 1 Sensor box
- 1 MiniDIN 8-pin connector cable, 60 cm long
- 1 Silicone hose, internal diameter 2 mm, 1 m long

#### 4. Technical data

Measurement range:	$\pm 1000$ hPa
Sensor type:	Semiconductor sensor
Accuracy:	$\pm 1$ %
Resolution:	$\pm 100$ hPa
Connections:	2 hose connections, 4.8 mm diameter

#### 5. Instructions

- Cut the silicone hose into sections of the required length.
- Using the lengths of hose, make the pressure connections between the sensor box and the Stirling engine.

- Note the “positive” and “negative” labelling of the hose connections - connect the hoses correctly according to the effective direction of the pressure.
- During the experiment, check that no elastic expansion of the hose is occurring – this can cause the pressure reading to be lower than the correct value.

## 6. Application

Measurement of the pressure difference in the transparent Stirling engine U10050, and analysis of the data using 3B NETlab™.

## 7. Sample experiment

### Recording operating pressures in Stirling engine U10050 while it is in motion

Apparatus required:

1 3B NETlog™ unit	U11300
1 3B NETlab™ program	U11310
1 Relative pressure sensor, $\pm 1000$ hPa	U11323
1 Transparent Stirling engine	U10050

- Set up the experiment as shown in fig. 1.
- Connect the relative pressure sensor to the 3B NETlog™ unit and wait for the sensor to be detected.
- Using a suitable length of silicone hose, make a pressure connection between the “positive” hose connection of the sensor box and the hose connection of the Stirling engine.
- Allow the engine to heat up and, after a few minutes, set it running.
- Open the application program (template) for the experiment with the  $\pm 1000$  hPa relative pressure sensor on the 3B NETlab™ unit.
- Measure the pressures.
- Evaluate the curve resulting from the measurements (fig. 2).


Fig. 1 Experiment set-up for recording operating pressures in Stirling engine U10050 while in motion


Fig. 2 Trace of pressure in Stirling engine U10050

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