Beschleunigungssensor 5-g U11363

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07/07 Hh

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07/07 Hh

1. Safety instructions

- Never exceed the maximum acceleration of 1500 times gravitational acceleration in any direction, to avoid permanent damage to the semiconductor built into the small black box!
- The maximum height from which the sensor can survive dropping onto a hard surface is 1.2 m.
- Only use the accelerometer 5-g for educational purposes!

The 5g accelerometer is not suitable for safety-related applications.

2. Description

Sensor box with permanently connected semiconductor acceleration sensor, Z-axis sensitive, for the measurement of gravity and the general acceleration of masses up to $\pm 5 \times g$.

The effective axis (Z-axis) is marked with an arrow and the label "Earth's Gravity Field" on the acceleration sensor.

"Capacitive" method of measurement (g-cell) with built-in signal linearisation, low-pass filtering, temperature compensation and automatic self-test.

The sensor box is automatically detected by the 3B NETlab™ interface.

3. Scope of delivery

- 1 Sensor box with permanently connected acceleration sensor, cable length 2 m.
- 1 Velcro strip, 500 mm long, 20 mm wide, self-adhesive
- 1 8-pin miniDIN connecting lead, 60 cm length
- 1 Instruction sheet for U11363

4. Technical data

Measurement range:	0 to ± 50 m/s ²
Sensor type:	Capacitive semiconductor sensor
Sensitivity:	Typically 400mV/g
Non-linearity:	No more than $\pm 1\%$ of the full measurement range.
Resolution:	0.03 m/s ²
Band width:	typically 50 Hz
Drill hole for sensor attachment:	3 mm diam. max.

5. Operation

- Place the sensor box alongside the experiment and attach the acceleration sensor (small black box) to the mass to be investigated (target). Use the supplied Velcro strip or a clamp for this purpose.
- Read the value of the acceleration from the display on the 3B NETlog™ unit.

6. Applications

Demonstration track and air track experiments:

- Downward acceleration
- Elastic and non-elastic impact

Oscillating spring-mass system

High-resolution measurement of objects' inclination

Pendulum oscillations

Jumping experiments; "bungee jumping"

7. Sample experiment

Acceleration measurement in a damped oscillating spring-mass system

Required equipment:

1 3B NETlog™ interface	U11300
1 3B NETlab™ software program	U11310
1 5-g accelerometer	U11363
1 Stand base	U13270
1 Stand rod, 750 mm length	U15003
1 Stand rod, 250 mm length	U15001
2 Universal clamps	U13255
1 Coil spring 3 N/m	U15027
1 Weight 100 g, from	U30016

- Set up the equipment for the experiment as in Fig. 1.

- Run the 3B NETlab™ software with the appropriate template for the experiment using the 5g accelerometer.
- Attach the acceleration sensor to the weight with a piece of Velcro.
- Suspend the weight and acceleration sensor from the eye at the bottom of the coil spring and be careful not to hinder the oscillating motion.
- Drape the connecting lead for the acceleration sensor over the universal clamp, as shown in Fig. 1. This adds further to the damping.
- Pull down the weight by hand to the level of the stand base and release it.
- Start recording the measurement data in 3B NETlab™ (Fig. 2).
- Analyse the recorded chart.

Fig. 1 Acceleration measurement for a damped oscillation of a mass on a spring


Fig. 2 Monitor display of the damped oscillation of a mass on a spring in 3B NETlab™ (U11310)

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