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
12/12 ALF

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### 1. Safety instructions

Hot cathode tubes are thin-walled, highly evacuated glass tubes. Treat them carefully as there is a risk of implosion.

- Do not subject the tube to mechanical stresses.
- Do not subject the cable connected with the anode to any tension.
- The tube may only be used with tube holder S (1014525).

If voltage or current is too high or the cathode is at the wrong temperature, it can lead to the tube becoming destroyed.

- Do not exceed the stated operating parameters.
- When the tube is in operation, the terminals of the tube may be at high voltages with which it is dangerous to come into contact.
- Only use safety experiment leads for connecting circuits.
  - Only change circuits with power supply equipment switched off.
  - Only exchange tubes with power supply equipment switched off.

When the tube is in operation, the stock of the tube may get hot.

- Allow the tube to cool before dismantling.

The EC directive on electromagnetic compatibility is only guaranteed when using the recommended power supplies.

### 2. Description

The dual beam tube can be used to determine the specific charge  $e/m$  from the diameter of the path followed by electrons fired into the tube from a perpendicularly mounted gun with a vertically aligned magnetic field and observation of the spiral path followed by electrons fired axially into a co-axial magnetic field.

The dual beam tube is a partly evacuated electron tube, filled with helium at low pressure and equipped with both axial and perpendicular electron guns. The electron beams are perpendicular to one another and a common deflector plate is provided for both guns. The electron beam source is an oxide cathode heated indirectly via a heating coil. The electron paths show up as a fine, slightly greenish beam due to impact excitation of the helium atoms.

### 3. Technical data

Filament voltage:	7.5 V AC/DC max.
Anode voltage:	100 V DC max.
Anode current:	30 mA max.
Deflector voltage:	50 V DC max
Glass bulb:	130 mm dia. approx.
Total length:	260 mm approx.
Gas filling:	Helium at 0.1 torr pressure

### 4. Operation

To perform experiments using the dual beam tube, the following equipment is also required:

1 Tube holder S	1014525
1 Power supply 500 V (115 V, 50/60 Hz)	1003307
or	
1 Power supply 500 V (230 V, 50/60 Hz)	1003308
1 Helmholtz pair of coils S	1000611
1 Analogue multimeter AM50	1003073

#### 4.1 Setting up the tube in the tube holder

The tube should not be mounted or removed unless all power supplies are disconnected.

- Press tube gently into the stock of the holder and push until the pins are fully inserted. Take note of the unique position of the guide pin.

#### 4.2 Removing the tube from the tube holder

- To remove the tube, apply pressure with the middle finger on the guide pin and the thumb on the tail-stock until the pins loosen, then pull out the tube.

### 5. Example experiments

#### 5.1 Determination of $e/m$

An electron of charge  $e$  moving at velocity  $v$  perpendicularly through a magnetic field  $B$  experiences a force  $F$  that is perpendicular to both  $B$  and  $v$  and the magnitude of which is given by:

$$F = evB$$

This causes the electron to follow a circular electron path in a plane perpendicular to  $B$ . The centripetal force for an electron of mass  $m$  is

$$F = \frac{mv^2}{R} = evB$$

which implies

$$B = \frac{v}{\frac{e}{m}R} \text{ tesla}$$

Rearranging the equation gives

$$\frac{e}{m} = \frac{v}{BR}$$

If the beam is subjected to a known magnetic field of magnitude  $B$ , and  $v$  and  $R$  are both calculated then the ratio  $e/m$  can be determined.

The law of conservation of energy means that the change in kinetic energy plus the change in potential energy of a charge moving from point 1 to point 2 is equal to zero since no work is performed by external forces.

$$\left(\frac{1}{2}mv_2^2 - \frac{1}{2}mv_1^2\right) + (eU_2 - eU_1) = 0$$

The energy of an electron in the dual beam tube is given by:

$$eU_A = \frac{1}{2}mv^2$$

By solving for  $v$  and replacing it in the equation

$$\frac{e}{m} = \frac{v}{BR}$$

the following emerges

$$\frac{e}{m} = \frac{2U_A}{B^2R^2}$$

The term  $e/m$  is the specific charge of an electron and has the constant value  $(1.75888 \pm 0.0004) \times 10^{11}$  C/kg.

#### 5.1.1 Determination of $B$

The Helmholtz coils have a diameter of 138 mm and give rise to a magnetic flux in Helmholtz configuration as given by

$$B = \mu_0 H = (4.17 \times 10^{-3}) I_H \text{ tesla}$$

and

$$B^2 = 17.39 \cdot 10^{-6} I_H^2$$

where  $I_H$  is the current in the Helmholtz coils.

The following are also true

$$\frac{e}{m} = \frac{U_A}{I_H^2 R^2} \cdot 1.15 \cdot 10^5 \text{ and}$$

$$I_H^2 = k \frac{U_A}{R^2}$$

#### 5.1.2 Determination of $R$

Referring to the diagram Fig. 1, the beam emerges from the electron gun at C travelling along the axis of the tube. The electron is then deflected in a circular path with the tube axis forming a tangent. The centre of this circle is at B and it lies in the plane of DCD' about 2 mm behind the plane of EE'.

$$AB^2 = BC^2 + AC^2 - 2BC \cdot DC$$

$$R = BC = AB = \frac{AC^2}{2DC} = \frac{x^2 + y^2}{2y}$$

$$R^2 = \left[ \frac{x^2 + y^2}{2y} \right]^2$$


Fig. 1 Derivation of  $R$

- Connect up the tube as in Fig. 3.
- Dim the room lighting.
- Assure the correct positioning of the Helmholtz coils in Helmholtz configuration.
- Set the heater voltage  $U_F$  to 7 V and wait about 1 minute for the heater temperature to achieve thermal stability (see remarks in section 7).
- Set the anode voltage  $U_A$  to 90 V (plate voltage  $U_P = 0$  V).
- Set the current in the coils  $I_H$  so that the deviated beam passes through point A on the edge of the fluorescent screen of the tube. Simultaneously focus the beam using a plate voltage  $U_P$  of no more than 6 V.
- Mark point A on the tube using a felt-tip pen.
- Increase  $U_A$  and set  $I_H$  so that the deflected beam always passes through A. Enter all the values into a table.

$U_A$ in volts	$I_H$ in amps	$I_H^2$
90		
100		
110		
120		

- Increase  $I_H$  so that the deflected beam always passes through point E and enter the values in a corresponding table.
- Mark point E on the tube using a felt-tip pen.
- Plot the graphs of the values from both tables

- Use a vernier calliper to measure the diameters  $AA'$ ,  $EE'$  and distance  $AE$ .
- Complete the table and calculate  $R^2$ .

AE mm	x = AE/ 2 mm	$x^2$ mm <sup>2</sup>	2y = EE' mm	y = EE' / 2 mm	$y^2$ mm <sup>2</sup>	$R^2$ mm <sup>2</sup>

- Replace the values in the equation

$$\frac{e}{m} = \frac{U_A}{I_H^2 R^2} \cdot 1.15 \cdot 10^5$$

and calculate a mean value for  $e/m$ .

## 5.2 Deflection in a circular path and the determination of $e/m$

- Connect up the tube as in Fig. 4.
- Assure the correct positioning of the Helmholtz coils in Helmholtz configuration.
- Set the anode voltage  $U_A$  to 100 V (plate voltage  $U_P = 0$  V).
- Set the current in the coils  $I_H$  so that the deflected beam moves in a circular path with the plane  $AA'$  tangential to it.

It is practical in this instance to observe the beam from above, from where it appears as a straight line and can be focused using  $U_P$  to a maximum of 6 V.

Note: the axial non-linearity of the beam has the effect of pushing the beam out of the plane of the electron gun. In order to obtain more accurate results, the tube should be turned within the brace that holds it so that the circular path is in the plane of the gun.  $I_H$  should also be modified so that plane  $AA'$  makes a good tangent with the path. A slight angle to the axis of the tube is tolerable. The beam travels in a slightly spiraling path instead of an accurate circle.

- Increase  $U_A$  and set  $I_H$  so that the plane  $AA'$  always forms a tangent to the deflected beam. Tabulate  $I_H$  against  $U_A$  and plot the graph.
- Evaluate  $R = AE/2$  and  $R^2 = AE^2/4$  like in experiment 5.1.
- Replace the values in the equation

$$\frac{e}{m} = \frac{U_A}{I_H^2 R^2} \cdot 1.15 \cdot 10^5$$

and calculate a mean value for  $e/m$ .

### 5.3 The effect of an axial magnetic field

- Connect up the tube as in Fig. 5.
- Place one coil into the groove of the tube holder from the front so that it encircles the screen and connect it to the 12 V output of the power supply.
- Set the anode voltage  $U_A$  to no more than 60 V (plate voltage  $U_P = 0$  V).
- Gradually increase the coil current  $I_H$ .

With only one axial velocity vector  $v_a$  the axial non-linearity of the beam is corrected and coincides with the true axis of the field.

- Mark the position of the beam with a felt-tip pen.
- Set  $I_H$  to 1.5 A and increase  $U_P$  gradually so that a second velocity vector  $v_p$  affects the beam.
- View the beam end-on through the coil.

The beam path turns into a helix. The beam no longer goes around the axis of the field but returns to a different position along the axis after every loop.


Fig. 2 Helical path of the deflected beam

- Reverse the polarity of the magnetic field  $B$  by reversing that of the Helmholtz coils and observe what happens to the beam.
- Change the anode voltage and observe its effect on the helical path of the beam. Then restore the voltage to 60 V.

### 6. Errors in the results

1. The circular beam path in experiment 5.2 is visible because of photo-emission. The energy involved in this process is lost and never replaced. This means that the beam actually tends to follow a spiral path instead of a circle. For a fixed radius  $R$  and a real circle  $U_A/I_H^2$  would be larger than the values that we measure. For this reason the error in the value of  $e/m$  is always on the negative side. Nevertheless results can be achieved that are accurate to within 20%.

2. In experiments where the beams are deflected into semi-circular paths as in experiment 5.1 results are larger than the published value. Points A and E, through which the beam is deflected, lie outside the homogeneous region of the Helmholtz coils so that the magnetic flux is reduced at these points. For a fixed radius  $R$  and a truly homogeneous field  $U_A/I_H^2$  would be smaller than the values we measure. For this reason the error in the value of  $e/m$  is always on the positive side. Nevertheless results can be achieved that are accurate to within 20%.

### 7. Remarks

1. Limiting of anode current: to avoid an excessive degree of emission of positive ions towards the electron emitting chemicals of the cathode, the anode current should be limited to below 20 mA wherever possible. Higher current may be tolerated for a short time but over long periods it reduces the lifespan of the tube.
2. Thermal stability of the cathode: for the same reason, you should avoid starting the electron gun when the cathode is cold and only just heating up.
3. Focussing the beam: Small voltages  $U_P$  applied to the deflector plates enable the beam to be focussed. Voltages greater than 6 V cause results to deteriorate.


Fig. 3 Determining  $e/m$  using the axial electron gun


Fig. 4 Determining  $e/m$  using the perpendicular electron gun



Fig. 5 The effect of an axial magnetic field

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