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07/12 ALF

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### 1. Safety instructions

Caution. Lamps get hot when they are switched on for a long period of time.

- When handling the lamp during the experiment, use a cloth or other protection. Do not cover the air slits.
- After the experiment leave the lamp to cool.

### 2. Description

The five-beam optical halogen lamp is a bright light source with five parallel light apertures for experiments involving ray optics to be conducted on a lab bench.

It consists of a metal housing with integrated ventilation fan and includes an adjustable reflecting mirror for setting the beam length. Thanks to a magnetic foil the lamp can also be used on a whiteboard.

### 3. Technical data

|                     |                                    |
|---------------------|------------------------------------|
| Halogen lamp:       | 12 V, 50 W                         |
| Connection:         | via 4-mm safety jacks              |
| Slit width:         | 2 mm                               |
| Slit spacing:       | 18 mm                              |
| Housing dimensions: | approx. 210x118x85 mm <sup>3</sup> |

### 4. Operation

#### 4.1 General notes

Suitable voltage sources for operation of the halogen lamp are, for example:

Transformer 12 V, 60 W (230 V, 50/60 Hz)  
1000593  
or  
Transformer 12 V, 60 W (115 V, 50/60 Hz)  
1006780

- Do not allow the lamp to suffer mechanical shocks.
- Do not supply the lamp with an operating voltage in excess of 12 V.

## **4.2 Accessories**

|                               |         |
|-------------------------------|---------|
| Set of optical components     | 1002993 |
| Spare halogen lamp 12 V, 50 W | 1002837 |

## **4.3 Changing the bulb**

Any deposits of fat from the skin on a halogen lamp bulb cause the glass to become fogged and significantly reduce the life of the bulb.

- Do not touch the glass bulb of the halogen lamp with fingers.
- Remove the securing screws and remove the lid of the casing.
- Take out the defective halogen bulb.
- Holding the replacement bulb with a cloth or other protective material, push the connecting wires into the socket.
- Screw the lid back on.

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