

# EMP 2 PRO

Dual-channel transcutaneous nerve stimulator



Art. No. 104026



Operating Instruction Art. No. 100846

**Content**

Purpose .....	3
Important safety instructions .....	3
Explanation of symbols .....	4
Device description .....	4
<b>Operating the EMP 2 PRO</b>	
Switching on the device .....	5
Activating the editing mode .....	5
Resetting the operating parameters and deleting the user programmes .....	5
Programme selection .....	5
Starting stimulation .....	5
Setting the intensity .....	6
Symbols on the display .....	6
Pause the stimulation .....	6
Stopping stimulation .....	6
Locking the programme .....	7
Locking the intensity .....	7
Activating the editing mode .....	7
Editing standard or user programmes .....	7
Programme parameters .....	8
Deleting the user programmes .....	9
Description of parameters .....	9
Calling up the parameters .....	9
Resetting the operating parameters .....	9
Switching off the device .....	9
<b>Description of the programmes .....</b>	<b>10</b>
<b>General information</b>	
Specifications .....	17
Changing the battery .....	18
Readjustments, alterations and repairs .....	18
Circuits diagrams, spare parts list and setting instructions .....	18
Warranty .....	18
Cleaning and care of the device .....	18
Classification .....	18
Technical checks .....	19
Combination .....	19
<b>Accessories</b>	
Rubber electrodes .....	19
Self-adhering electrodes .....	21
Accessories .....	22
<b>Electrode position</b>	
programmes pain therapy P2-P4-P5-P6-P7-P8-P10 .....	23
programmes relaxation P1, P3, P9 and P20 .....	24
programmes for muscle stimulation and re-education P11-P16 .....	25

**Purpose**

The EMP 2 PRO device is designed for transcutaneous muscle and nerve electrostimulation. Please do not use this device for any other purpose.

**Important safety instructions**

Please read the „Operating Instruction“ carefully before using the device.

**To ensure safe use of the device**

- Only use the EMP 2 PRO with original accessories.
- Keep the EMP 2 PRO away from water or other fluids.
- Do not drop the EMP 2 PRO, handle it incorrectly or expose it to extreme temperatures or high humidity (only use at temperatures between 10 °C and 40 °C and at a relative humidity below 90 %).
- Never use the EMP 2 PRO when it is not functioning properly or when it is damaged in any way.
- Store the EMP 2 PRO in its original packaging after use to protect it from damage and contamination.

**Precautions**

Patients with an implanted electronic device (e. g. pacemaker) should not undergo electrotherapy with the EMP 2 PRO before consulting a doctor. The EMP 2 PRO may only be connected to one person at a time.

**Warning**

Connecting the patient at the same time to a high-frequency surgical unit may cause burns under the electrodes.

Operation near (e. g. 1 m) from a short-wave, mobile phones or microwave device may cause fluctuations in the baseline values of the electrotherapy device.

The minimum area of the electrodes should not be less than 2 cm<sup>2</sup>.

## Explanation of symbols



Attention: Read Operation Instructions!



BF type application part. Protection against electric shock.



The year of construction of the product follows this symbol.

**REF**

The article or order number of the product follows this symbol.

**SN**

The serial number of the product follows this symbol.



This equipment is marked with the recycling symbol. It means that at the end of the life of the equipment you must dispose of it separately at an appropriate collection point and not place it in the normal domestic unsorted waste stream. This will benefit the environment for all.

CE 0197

Conform with Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

## Device description

The EMP 2 PRO was designed to stimulate human nerves and muscles.

The keys may be used to make all adjustments necessary. The display shows the different operating modes.

1. Display
2. Programme key
3. Key for visualisation of the parameters, pause and edit mode
4. Modification keys (these keys may be used to alter intensity during stimulation and to alter parameters in the editing mode)
5. ON/OFF key
6. Battery compartment
7. OUT hubs
8. 230-V transformation socket



## Operating the EMP 2 PRO

### Switching on the device

Switch on the device using the key. The EMP 2 PRO starts with the programme number which was active when the device was last switched off.


### Using the 230-Volt mains supply (optional)

After having connected the transformer to your 230-Volt socket, connect the cable to the central supply socket in the front face of the appliance between the channel outlets 1 and 2. It should be noted that when the machine is connected to the transformer connected to the 230-Volt mains supply it will operate from the mains supply and not the battery.

### Activating the editing mode


Only if you switch on the EMP 2 PRO pressing the **E** key, you can activate the standard programmes (please refer to section entitled "Activating the editing mode") and make copies of them in the form of user programmes. However, this function is only available if the device was not locked beforehand.

### Resetting the operating parameters and deleting the user programmes



If you wish to delete the user programmes and return the operating time meter, the cycle meter and the mean stimulation intensity to zero, keep the **▼** key on the left and the **P** key pressed down while switching on the device (please refer to section entitled "Resetting the operating parameters").

### Programme selection



Use the **P** key to change to the next standard or user programme. However, this is only possible with an unlocked device. The standard programmes 1 to 20 marked with a "P" run initially, followed by the user programmes marked with an "U", if present. When the last programme has been reached, pressing the key again returns the setting to programme 1 (P1).

### Starting stimulation

First, position the electrodes at the desired points of the body and connect the electrode cables to the electrodes and the device. When the desired programme has been selected using the **▲▼** modification key, or the **E** key increase the intensity (Fig. 2), the two keys may be used to start stimulation. In





the display on the right hand-side you find two vertical lines. This means: with the **E** key both channels will be increased at the same time.

### Setting the intensity

You can use the **▲ ▼** modification keys or the **E** key to set the intensity of both channels to a pleasant value at any time. The display shows the strength of current for both channels, which may be varied from 0-100 mA. Intermediate values are clearly marked by a dash behind the digit. If the electrodes are not correctly connected to the device, the intensity is returned to zero starting from a current of 4 mA. If you are in the run-down or pause phase in muscle stimulation programmes and you press one of the modification keys, the working phase is started directly from the build-up phase.

### Symbols on the display

During stimulation, the current intensity status is shown in the middle of the display by four symbols:

-  Continuous stimulation or work phase in muscle stimulation. Stimulation is performed at the set intensity.
-  Build-up phase in muscle stimulation. Intensity is brought up to the set value from zero. In dynamic stimulation, the intensity of channel 1 is adjusted upwards to the set value, whereas that of channel 2 is adjusted downwards to zero.
-  Run-down phase in muscle stimulation. The intensity is brought down from the set value to zero. In dynamic stimulation, the intensity of channel 2 is adjusted upwards to the set value, whereas that of channel 1 is adjusted downwards to zero.
-  Pause phase in muscle stimulation. Stimulation is performed with 75 % of the set intensity at a frequency of 3 Hz.

### Pause the stimulation




When the sequence number is not flashing, the **E** key may be used to interrupt stimulation. The word "PAUSE" appears on the display. If the pause key is pressed again, the last intensity set for both channels is restored via a build-up ramp.

### Stopping stimulation

Use the **P** key or modulation keys to stop stimulation at any time. Stimulation stops automatically at the end of a programme.

### Locking the programme




Choose the programme desired for the patient using the **P** key (please refer to section entitled "Description of programmes"). Lock the device by then pressing the **▼** key on the left and also the **P** key. A key appears in the upper right-hand corner of the display. All device functions apart from stimulation with the selected programme are now no longer available. Use the same key combination to unlock the device. A locked device cannot be started in the editing mode.

### Locking the intensity

Choose the programme desired for the patient using the key (please refer to section entitled "Description of the programmes"). The EMP 2 PRO will be locked automatically after 10 sec. without increasing or decreasing the intensity. On the display in the right corner the "P" symbol will be showed. To unlock the EMP 2 PRO you have to decrease the intensity.

### Activating the editing mode



Please switch on the device while pressing the **E** key. An "E" on the right-hand side of the display indicates the editing mode. Then use the key to select the standard or user programme (Fig. 6) that you wish to edit.

The key is then used to start editing the currently active standard or user programme. This is shown by the fact that the "P" or "U" symbol before the programme number changes to an "E".

### Editing standard or user programmes

The parameters of each programme are used as the starting value. Press the **E** key to jump to the next parameter or to save the values.

Use the **P** key to interrupt the editing process at any time without saving the altered parameters.

Various parameters may be processed, depending on the type of programme. If not indicated otherwise, the frequency and the pulse width apply for both channels. Use the modification keys to adjust the values. If you keep the keys pressed down, counting up or down of the parameter is automatically continued.

## Programme parameters



### Programme 1, 3, 4, 8, 9, 10, 11:


1. frequency: 20 Hz up to 120 Hz (Fig. 7)
2. pulse width: 70  $\mu$ s up to 500  $\mu$ s (Fig. 8)
3. total run time: 1 up to 99 min. (Fig. 12)



### Programme 2, 5 (TENS):

#### channel 1 gate control + channel 2 endorphine release

1. frequency channel 1: 80 Hz, 90 Hz, 100 Hz, 110 Hz, 120 Hz (Fig. 9)
2. frequency channel 2: 0,5 Hz, 1 Hz, 2 Hz, 5 Hz, 10 Hz
3. pulse width: 60  $\mu$ s up to 500  $\mu$ s (Fig. 8)
4. total run time: 2 up to 99 min. (Fig. 12)




### Programme 6:

1. total run time: 1 up to 99 min. (Fig. 12)

### Programme 7 :

#### Epicondylites modulation


1. max. frequency: 1 Hz up to 120 Hz (Fig. 10)
2. min. frequency: 0,5 Hz up to max. frequency -8 (Fig. 11)
3. pulse width: calculated according to the frequency
4. total run time: 1 up to 99 min. (Fig. 12)



### Programme 12, 13, 14, 15, 16, 18, 19:

#### muscle stimulation


1. frequency: 20 Hz up to 120 Hz (Fig. 7)
2. pulse width: 70  $\mu$ s up to 500  $\mu$ s (Fig. 8)
3. Ramp up: 1 to 30 s
4. Ramp down: 1 to 30 s
5. Work time: 1 to 30 s
6. Total run time: 1 up to 99 min. (Fig. 12)



### Programme 21:

#### muscle stimulation

1. frequency: 20 Hz up to 120 Hz (Fig. 7)
2. pulse width: 60  $\mu$ s up to 500  $\mu$ s (Fig. 8)
3. Ramp up: 0,5 to 4 s
4. Ramp down: 0,5 to 2 s
5. Work time: 1 to 25 s
6. Total run time: 1 up to 99 min. (Fig. 12)




To save the new programme, press the **E** key once more after processing the last parameter. An "U" symbol now appears before the programme number on the display.



## Deleting the user programmes

To delete all user programmes, keep the **P** key and the **▼** key on the left pressed down when switching on the device. This will also return all operating parameters to zero.

## Description of parameters



These values are used to check the patient's stimulation behaviour. The following information are reported by the EMP 2 PRO: The stimulation time (Fig. 13) is the total time during which stimulation was carried out with the device. The cycles (Fig. 14) mark how often the device was switched on and concomitantly used for stimulation. The mean stimulation intensity (Fig. 15) is the mean value of all values for current set by the patient. This information is registered separately for the two channels. The values are updated every minute during stimulation. However, only intensities above 4 mA are considered.



## Calling up the parameters

Press the **▼** key on the left and the **E** key at the same time. The stimulation time in hours and minutes is first shown in the upper part of the display. The number of cycles is shown in the lower part of the display. After pressing the **E** key, the mean stimulation intensity for channel 1 is shown in the lower left part of the display, that for channel 2 in the lower right part.

## Resetting the operating parameters

To reset all the described values to zero, keep the **▼** key on the left and the **P** key pressed down when switching on the device. Please be aware that all user programmes are also deleted.

## Switching off the device

Use the **⊙** key to switch off the device. If the battery voltage has fallen below a critical value, or if no key is pressed for two minutes outside stimulation, the device switches itself off automatically. This is signalled by a beep.

## Description of the programmes

### Programme 1 Lumbalgie gate control


**Type of stimulation:** Gate Control 80 Hz and dynamic stimulation

**Parameters:** Frequency channel 1 and 2 = 80 Hz  
Pulse width = 150  $\mu$ s  
Treatment time = 20 min

**Symbol showing on display during stimulation:**

Ramp up: ▲ Ramp down: ▼

**Description:** Both channels are continuously operated with the same frequency and pulse width.




### Programme 2 Lombo-sciatalgie

**Type of stimulation:** Gate Control 80 Hz (channel 1)  
+ 2 Hz endorphin release (channel 2)

**Parameters:** Frequency channel 1 = 80 Hz  
Frequency channel 2 = 2 Hz  
Pulse width = 200  $\mu$ s  
Treatment time = 20 min

**Symbol showing on display during stimulation:** ■

**Description:**



### Programme 3 Cervicobrachial neuralgia


**Type of stimulation:** Gate control 80 Hz stimulation fast dynamic movement

**Parameters:** Frequency = 80 Hz  
Pulse width = 150  $\mu$ s  
Treatment time = 20 min  
Ramp up = 0,2 s  
Ramp down = 0,2 s

**Symbol showing on display during stimulation:**

Ramp up: ▲ Ramp down: ▼

**Description:** Basically, the device works like in the muscle mode but without working and break times. During the build-up phase intensity is increased up to the set value. The run-down phase is following immediately, decreasing intensity back to zero. The channels of the EMP 2 PRO work alternately, i. e. when channel 1 has reached its maximum intensity; the intensity of channel 2 is zero and vice versa. There is a convenient massage effect owing to the dynamic stimulation when placing the electrodes next to each other on the painful area.

**Programme 4 Chronic pain****Type of stimulation:** Endorphin release Burst 2 Hz**Parameters:** Frequency = 80 HzPulse width = 100  $\mu$ s


Treatment time = 30 min

Work time = 0,25 s

Pause = 0,25 s

**Symbol showing on display during stimulation:**

Working time: ■ Pause: ■

**Description:****Programme 5 Gate control + endorphin release****Parameters:** Frequency channel 1 = 100 Hz


Frequency channel 2 = 2 Hz

Pulse width = 200  $\mu$ s

Treatment time = 20 min

**Symbol showing on display during stimulation: ■**

**Description:** Channel 1 works continuously at a high frequency between 80 and 120 Hz, whereas channel 2 operates at a low frequency between 0,5 and 10 Hz.



## Programme 6 Gonarthrosis-Coxarthrosis

**Type of stimulation:** HAN


**Parameters Phase 1:** Frequency = 100 Hz  
Pulse width = 100  $\mu$ s  
Work time = 3 s

**Parameters Phase 2:** Frequency = 2 Hz  
Pulse width = 200  $\mu$ s  
Work time = 3 s  
Treatment time = 20 min

**Symbol showing on display during stimulation:**

Phase 1:  Phase 2:

**Description:** The device stimulates with a high frequency and a relatively low pulse width (phase 1) respectively with a low frequency and a relatively high pulse width (phase 2) alternately. The duration of the phases is 3 s each.



## Programme 7 Epicondylites


**Type of stimulation:** Modulation

**Parameters:** Frequency minimal = 2 Hz  $\Rightarrow$  Pulse width = 200  $\mu$ s  
Frequency maximal = 80 Hz  $\Rightarrow$  Pulse width = 100  $\mu$ s  
Treatment time during modulation = 7,5 s  
Treatment time total = 30 min

**Symbol showing on display during stimulation:**

Ramp up:  Ramp down:

**Description:** In the build-up phase, frequency is increased within a few seconds in individual steps, starting from minimum and reaching maximum frequency, the pulse width being adapted accordingly. This is followed by a run-down phase in which the intensity is regulated back down to the minimum.



## Programme 8 Algodystrophy


**Type of stimulation:** Modulation

**Parameters:** Frequency minimum = 2 Hz => Pulse width = 148  $\mu$ s  
 Frequency maximum = 80 Hz => Pulse width = 70  $\mu$ s  
 Treatment time during modulation = 7,5 s  
 Treatment time total = 20 min

**Symbol showing on display during stimulation:**

Ramp up:  Ramp down: 

**Description:** In the build-up phase, frequency is increased within a few seconds in individual steps, starting from minimum and reaching maximum frequency, the pulse width being adapted accordingly. This is followed by a run-down phase in which the intensity is regulated back down to the minimum.



## Programme 9 Gate control dynamic anti-adaptation


**Type of stimulation:** Gate control 80 Hz stimulation dynamic

**Parameters:** Frequency = 80 Hz  
 Pulse width = 150  $\mu$ s  
 Treatment time = 20 min  
 Ramp up = 0,5 s  
 Ramp down = 0,5 s

**Symbol showing on display during stimulation:**

Ramp up:  Ramp down: 


**Description:**



## Programme 10 Gate control classic 80 Hz

**Parameters:** Frequency = 80 Hz  
 Pulse width = 150  $\mu$ s  
 Treatment time = 20 min




**Symbol showing on display during stimulation:** 

**Description:****Programme 11 Treatment of contractures**

**Parameters:** Frequency = 1 Hz  
 Pulse width = 150  $\mu$ s  
 Treatment time = 20 min

**Symbol showing on display during stimulation:** See programme 10

**Muscle stimulation programmes (P12-P16 and P18-P19)**


Ramp up:  Ramp down:  Pause: 

**Programme 12**

**Type of stimulation:** Atrophy upper extremities

**Type of stimulation:** Muscle strengthening

**Parameters:** Frequency = 35 Hz  
 Pulse width = 150  $\mu$ s  
 Ramp up = 2 s  
 Working period = 3 s  
 Pause = 9 s

**Description:****Programme 13 Atrophy upper lower extremities**

**Type of stimulation:** Muscle strengthening

**Parameters:** Frequency = 35 Hz  
 Pulse width = 300  $\mu$ s  
 Ramp up = 2 s  
 Working period = 3 s  
 Pause = 9 s

**Programme 14 Muscle strengthening upper extremities****Type of stimulation:** Muscle strengthening**Parameters:** Frequency = 65 Hz  
Pulse width = 250  $\mu$ s  
Ramp up = 2 s  
Working period = 4 s  
Pause = 8 s

---

**Programme 15 Muscle strengthening lower extremities****Type of stimulation:** Muscle strengthening**Parameters:** Frequency = 65 Hz  
Pulse width = 300  $\mu$ s  
Ramp up = 2 s  
Working period = 4 s  
Pause = 8 s

---

**Programme 16 Venous reflux****Parameters:** Frequency = 35 Hz  
Pulse width = 250  $\mu$ s  
Ramp up = 3 s  
Working period = 5 s  
Pause = 10 s

---

**Programme 17 URO vesical instability****Parameters:** Frequency = 10 Hz  
Pulse width = 180  $\mu$ s  
Treatment time = 20 min**Description:** See Programme 10

---

**Programme 18 URO incontinence mixed****Parameters:** Frequency = 20 Hz  
Pulse width = 180  $\mu$ s  
Ramp up = 2 s  
Working period = 4 s  
Pause = 4 s

---


**Programme 19 URO incontinence urge****Parameters:** Frequency = 50 Hz  
Pulse width = 180  $\mu$ s  
Ramp up = 2 s  
Working period = 3 s  
Pause = 6 s

**Programme 20 Relaxation and stimulation dynamic****Type of stimulation:** Stimulation dynamic

**Parameters:** Frequency = 120 Hz  
 Pulse width = 150  $\mu$ s  
 Treatment time = 20 min  
 Ramp up = 2 s  
 Ramp down = 2 s

**Symbol showing on display during stimulation:**

Ramp up:  Ramp down: 

**Description:****Programme 21 Agonist/Antagonist****Type of stimulation:** Stimulation dynamic - channel 1, channel 2

**Parameters:** Frequency = 40 Hz  
 Pulse width = 300  $\mu$ s  
 Treatment time = 20 min  
 Ramp up = 2 s  
 Ramp down = 1 s  
 Working period = 3 s  
 Pause = 6 s

**Description:** First increase intensity channel 1 when the contraction is there you can increase channel 2 and when there is a contraction the treatment will start. When you want to increase channel 1 during stimulation channel 2 the intensity of channel 2 will go to zero and you can increase channel 1.



## General information


### Specifications

Dual-channel nerve stimulator with electrically insulated channels, constant current characteristic, output short circuit control element (AKS) and 20 changeable integrated programmes.


Output current	100 mA (with 1 k $\Omega$ real)
Frequencies	0,5-120 Hz
Pulse width	70-500 $\mu$ s
Current consumption	12 mA
Current supply	3 x 1,5 V compound battery
Dimensions approx.	132 mm x 83 mm x 39 mm
Weight approx.	190 g (without accessories)

### Pulse forms


**Pulse form:  
with 1k $\Omega$  load real**



**With ANSI/AAMI  
standard load**



**Change in output current  
depending on load**



## Changing the battery

The voltage of the batteries in the device is displayed during operation by the segments within the battery symbol. If the voltage falls below a critical value, the EMP 2 PRO switches itself off automatically and cannot be switched back on again. You then need to insert new batteries into the device.

- Switch off the device
- Take off the lid of the battery compartment
- Remove the used batteries from the compartment
- Insert new batteries into the compartment

### **PLEASE DISPOSE OF USED BATTERIES CORRECTLY.**

Please observe the instructions accompanying the battery charger if rechargeable batteries are used.

## Readjustments, alterations and repairs

The manufacturer is only responsible for the safety and performance of the EMP 2 PRO device when readjustments, alterations and repairs are carried out by authorised individuals and when the EMP 2 PRO is used in accordance with the user instructions.

## Circuits diagrams, spare parts list and setting instructions

Qualified technicians who are familiar with the technical features of the device can be provided with circuit diagrams, spare parts list and setting instructions by the manufacturer.

## Warranty

We give a guarantee of 1 year from the date of purchase on the EMP 2 PRO device. This guarantee does not cover cables and electrodes.

## Cleaning and care of the device

No special cleaning or care agents are required for the EMP 2 PRO. Clean the EMP 2 PRO with a soft, lint-free cloth. Please ensure that no moisture gets into the device. If moisture does enter the device, a technical check must be carried out before using the device again.

## Classification

In accordance with the Law on Medical Devices, the EMP 2 PRO device is classified as being a Class IIa medical device. Technical checks on the device should be performed every 24 months.

## Technical checks

These include:

1. Checking to see whether the user instructions are included in the accompanying documentation.
2. Checking the equipment for completeness.
3. Visual check:
  - a. for mechanical damage
  - b. for damage to all cables and connections
4. Functional safety
  - a. Checking the output signals with a load resistance of 1 k $\Omega$  real (current and voltage)
  - b. Checking the frequency
  - c. Checking the pulse width

These technical checks may only be performed by individuals with appropriate training. The results must be noted in the medical device book along with the date and name of the person carrying out the check.

## Combination

The EMP 2 PRO may be used together with all accessories mentioned in the chapter „Accessories“.

# Accessories

## Rubber electrodes

### Technical Data

Compound:	silicone-carbon
Durability:	approx. 12 months
Colour:	black
Manufacturer:	Pierenkemper GmbH, Germany










### Application

Spread 2-3 drops of schwa-medico electrode gel evenly over the flat side of the electrode. The electrode gel improves the conductivity of the electrode. Place the electrodes on the relevant skin areas and secure the electrodes with a tape strip. Do not apply to broken skin!

Because conductivity of the electrodes slowly decreases after approx. 6 months of use, replace electrodes after approx. 12 months of intensive use.

### Cleaning

Clean electrode surface with water and soap or appropriate disinfectant (e. g. alcohol 90 %) after every use.

Art. No.	Size	Quantity
107090	Rubber electrode, round 20 mm Ø 	2
107060	Rubber electrode, round 25 mm Ø 	2
107075	Rubber electrode, 40 x 28 mm 	2
107035	Rubber electrode, 56 x 28 mm 	2
107020	Rubber electrode, 75 x 30 mm 	2
107055	Rubber electrode, 90 x 35 mm 	2
107011	Rubber electrode, 38 x 45 mm 	2
107010	Rubber electrode, 48 x 48 mm 	2
107050	Rubber electrode, 70 x 65 mm 	2

Art. No.	Size	Quantity
107070	Rubber electrode, 70 x 140 mm	2



## Self-adhering electrodes

### Technical Data

Compound:	conductive and bonding material
Durability:	80-150 treatments
Colour:	grey
Manufacturer:	Pierenkemper GmbH, Germany

### Application

Place self-adhering electrodes directly on appropriate skin treatment area, no special preparation of electrode is necessary. See „Electrode maintenance“ below to insure electrode’s effectiveness. Do not apply to broken skin!

### Electrode maintenance





For hygienically reasons, self-adhering electrodes are for single patient use only. After every treatment session, return electrodes to their original foil and store them in their plastic bag. To increase longevity, store in cool area (i. e. refrigerator). With appropriate use and proper care the electrodes are applicable 80-150 times. To renew adhesion of electrodes, place 2-3 drops of water on the electrode’s adhesive surface and air dry a couple of spronds before placing on treatment area.

Art. No.	Size	Quantity
281000	Stimex, round 32 mm Ø	4







282000	Stimex, round 50 mm Ø	4
--------	-----------------------	---



Art. No.	Size	Quantity
283400	Stimex, 50 x 50 mm	4
		
283600	Stimex, 50 x 90 mm	2
		
283000	Stimex, 50 x 130 mm	2
		
283100	Stimex, 80 x 130 mm	2
		

## Accessories

106712	Crocodile clips	4
		
108000	Electrode gel	60 g
		
106711	Electrode cable type 5	1
		
109000	Tape (type silk)	3 rolls
		

## Electrode position

### Electrodes position programmes pain therapy P2-P4-P5-P6-P7-P8-P10

P2 Lumbago (channel 1 and position Para vertebrae and the lumbar side, channel 2 direct of the ischiadicus nerve)

P4 Chronic pain burst endorphin release

P5 Gate control (channel 1 at 100 Hz) + endorphin release (channel 2 at 2 Hz)

P6 Gonarthrosis-Cxarthrosis

P7 Epicondylitis


P8 Algodystrophie

P10 Gate control classic

Neck




Deltoideus




Shoulder




Knee



Elbow



Trapezoid



**Electrode placement for dynamic stimulation for pain therapy and relaxation****Electrode placement for relaxation P1, P3, P9 and P20**

How the electrodes should be positioned for the programmes


P1 Lumbago

P3 Cervical brachial neuralgia

P9 Massage pain

P20 Muscle relaxation


Latissimus




Trapezoid



Para vertebra



Calf



Quadriceps



Leg flexor





**Electrode placement for muscle stimulation and re-education programmes P11-P16**

How the electrodes should be positioned for the programmes re-education

P11 Treatment of contractures

P12 Atrophy upper extremities

P13 Atrophy lower extremities

P14 Muscle training upper extremities

P15 Muscle training lower extremities

P16 Dynamic venous reflux

Hand flexor



Hand extensor



Biceps



Triceps



Back



Quadriceps




Leg flexor (hamstrings)



Foot flexor



m. Tibialis



Calf (Gastrocnemius)





**Informations:**

schwa-medico GmbH  
Export Department  
Wetzlarer Strasse 41-43  
35630 Ehringshausen · Germany  
Phone: +49 6443 8333-113  
Fax: +49 6443 8333-119  
export-sm@schwa-medico.com  
www.schwa-medico.com



schwa-medico  
Medizinische Apparate Vertriebsgesellschaft mbH  
Wetzlarer Straße 41-43 · 35630 Ehringshausen · Deutschland  
Tel. 06443 8333-444 · Fax 06443 8333-450

Manufacturer:  
Pierenkemper GmbH  
Hörnshheimer Eck 19 · 35578 Wetzlar · Germany  
www.pieren-tech.de  
Art.-Nr. 451600-0251

CE 0197