

3B SCIENTIFIC ACUPUNCTURE

Visit us at 3bscientific.com

ACUPUNCTURE NEEDLES



VINCO® Detox Acupuncture Needles

- Most efficient package for detox protocols with FingerSaver™ feature
- Metal pipe handles, 10 needles per pack (100 per box)
- Cost effective for detox clinics



Needle Length	Diameter	Item No.	Price
0.25 in	0.22	W70004	\$42.00
0.5 in	0.22	W70005	\$42.00
0.25 in	0.2	W70006	\$42.00
0.5 in	0.2	W70007	\$42.00



Perfect for:

- + General acupuncture
- + Quick needle access
- + Detoxification treatment

A. Vinco®-Cluster - Acu Needle 500 box

- 5 needles bundled in a SkinSaver™ guide tube
- Extended medical sealing paper for easy peeling
- 5 blister packs per flat with FingerSaver™ feature
- 20 flats, 500 needles per box

Needle Length	Diameter	Item No.	Price
1.0 in	0.3	W70033	\$23.50
1.5 in	0.3	W70034	\$23.50
0.5 in	0.25	W70035	\$23.50
1.0 in	0.25	W70036	\$23.50

B. Vinco®-EZY - Acu Needle 200 box

- All needles are separately positioned in a needle dispenser
- 20 needles plus 2 SkinSaver™ guide tube per dispenser
- 10 pouches, 200 needles per box

Needle Length	Diameter	Item No.	Price
2.0 in	0.3	W70076	\$13.50
0.5 in	0.25	W70077	\$13.50
1.0 in	0.25	W70078	\$13.50
1.5 in	0.25	W70079	\$13.50

C. Vinco® Press Tack - Press Tack 100 box

The needle base is anchored between two layers of surgical adhesive tape. Tack will not accidentally fall into ear canals or be left behind when removing it. The base is isolated from the skin and will not cause irritation. US-made breathable surgical tape, suitable for sensitive skin.

Needle Length	Diameter	Item No.	Price
1.50 mm	0.22	W70090	\$5.00
1.60 mm	0.2	W70091	\$5.00
0.9 mm	0.18	W70092	\$5.00
1.10 mm	0.18	W70093	\$5.00

➤ More sizes available at 3bscientific.com



VINCO® Blister - Acu Needle 100box

- Individually packaged with SkinSaver™ guide tubes
- 4mm needle length and guide tube differential
- 5 needles per flat with FingerSaver™ feature
- 20 flats, 100 needles per box

Needle Length	Diameter	Item No.	Price
25mm	0.2	W70026	\$7.25
40mm	0.3	W70012	\$7.25
50mm	0.2	W70028	\$7.25
75mm	0.22	W70024	\$7.25
13mm	0.16	W70031	\$7.25
25mm	0.25	W70016	\$7.25
40mm	0.25	W70017	\$7.25
50mm	0.25	W70018	\$7.25

VINCO® Intradermal, Acu Needle 100box

- Sterile
- Vacuum melted high purity stainless steel body
- Individually blister packaged
- Disposable, 100 needles per box



Needle Length	Diameter	Item No.	Price
6mm	0.12	W70008	\$7.00
3mm	0.14	W70009	\$7.00
4mm	0.14	W70010	\$7.00

➤ This pocket pen models is made in Japan and has an attractive chrome finish with a removable pocket clip.



Pulsar Piezo Stimulator

Generates a moderate stimulation pulse by pushing the activator at the end of the unit.

W70099 \$45.00



Stainless Steel Straight Probe

Probes for locating and treating ear or body points. 6" (15cm) long octagonal shaped handle with a 0.1" (0.3cm) ball-shaped tip for inducing sensation and applying pressure to treatment points.

W70096 \$6.50



D&D Acupuncture Needles (1000 needles/box)

- LIFTed™ needle body that offers maximum patient comfort
- Bright fluorescent orange handles for easy locating
- Thicker and longer plastic handles for easy handling
- Efficient package for detox protocols
- Great for detox clinics
- 100 flats, 10 needles per flat

Needle Length	Diameter	Item No.	Price
0.25 in	0.22	W70000	\$42.00
0.5 in	0.22	W70001	\$42.00
0.25 in	0.2	W70002	\$42.00
0.5 in	0.2	W70003	\$42.00

MOXA, EAR SEEDS & PELLETS



Onden Smokelss Moxa
 Kamaya Ondan Smokeless Needle Warming Moxa. Includes 6 needle caps, 1 spoon, 120/box.
W70111 \$48.50

MAXIM Ultra Pure Jr. Moxa
 This 4.6" x 0.25" slender tightly packed New Ultra Pure Jr. Moxa Rolls contain a finer grade of moxa. Burns at approximately 36 mins per stick, 20 sticks per box.
W70108 \$3.75

MAXIM EZY-Lit Smokeless
 This 4.5"x 0.55" Hollow Center EZY-Lit Smokeless moxa sticks are the same composition and size as the familiar Smokeless Moxa Sticks (MX-01). Burns at approximately 40 minutes per stick. 5 sticks per box.
W70109 \$4.00

MAXIM Smokeless Moxa
 This 0.5" x 4.5" Smokeless moxa sticks are the same composition and size as the familiar smokeless Moxa Sticks (MX-01). Use candle to ignite the roll. Burns at approximately 1.8" Per hour. 5 Sticks per box.
W70112 \$4.00

No Gold paint on the stick

Moxa Holder Extinguisher
 Made of chrome-plated solid brass. Can hold a burning moxa roll in the notch or extinguish most moxa rolls in the extinguish chamber in less than 15 seconds.
W70114 \$7.00

Moxa Essence Plus Spray - 100ml
W70115 \$11.00

Moxa Essence Plus Dropper - 30 ml
W70116 \$7.00


➤ More seeds & pellets available online at 3bscientific.com



A. Helio® Ear Seeds, Latex free (100 per box)
 Selected Vaccaria seeds approximately 2mm size and processed with rice wine for enhanced medicinal effect.
W70101 \$4.50

B. Helio® Magnetic Ear Pellets (100 per box)
 Magnetic pellets with random magnetic pole orientation, 1.7mm in diameter with tan colored surgical tape. Best applied with use of DP-01 probe.
W70102 \$5.00

+
Pre-taped:
 + Simply peel & apply!



A. Body Acupuncture Chart

This thickly laminated anatomical chart is printed on premium glossy UV resistant paper and comes with 2 sided lamination and metal eyelets. It can be written on and wiped off with non permanent markers. Size 38.5 x 26.7 in.

Laminated	VR1820L	\$21.00
Unlaminated	VR1820UU	\$14.00

B. Ear Acupuncture Chart

Laminated	VR1821L	\$21.00
Unlaminated	VR1821UU	\$14.00

Acupuncture Model

This new acupuncture model sets the standard in didactics and aesthetics and has been developed in cooperation with experienced acupuncture specialists.

It contains 396 acupuncture points. All conduction paths are colored according to the 5 phases.

We have intentionally dispensed with Chinese labelling on the model. This makes the accurate identification of acupuncture points significantly easier.

Male	N30	\$202.00
Female	N31	\$202.00



➤ 27.5 inches tall

SKINlike™ Ear set, one left and right ear

Practice without pain to perfect your technique. These original-sized replicas will facilitate learning and practicing ear acupuncture (also auricular medicine or auricular therapy).

The models are made of high-quality 3B SKINlike™ silicone – guaranteeing a long life and realistic ear surface structure. To aid in finding the acupuncture points, 18 anatomical structures are cited in an enclosed product manual. Supplied with 5 acupuncture needles.

3.7 x 2.4 x 1.6 in

N15 \$63.00

SKINlike™ Individual Acupuncture Ears

The same high quality replicas as item N15.

Acupuncture Ear - Left
N15/1L \$36.00

Acupuncture Ear - Right
N15/1R \$36.00

SKINlike™ Ears, set for 10 students

Includes:

- 10 Pairs of acupuncture ears (20 pcs)
- 50 Acupuncture needles
- 3.7 x 2.4 x 1.6 in

N16 \$600.00



3B LASER



Noninvasive, gentle and pain free 3B LASER

The innovative **3B LASER** is classified as **low-level laser therapy**, and offers **painless** and **noninvasive** treatment. Light particles (photons) with high excitation density are conducted through fine optical fibers. The fibers do not puncture the skin, but are placed on therapy points, radiating into the tissue. The lasers can be used simultaneously on any point combinations on the body, head, and ear. Thus, corresponding methodically to the **traditions of classical acupuncture**. They are equivalent to acupuncture treatments in their cerebral and peripheral stimulation effects.

Monitor the **3B LASER** via the mobile touch screen (**Control Touch Pad**), connected with the laser unit by a radio interface. You achieve optimum communication with your patients, who can follow the treatment at all times and thus be more actively involved in the elucidation and therapy. The daily work is greatly simplified so that both a bedside and remote observation of the therapy course is possible.

Optimised functionality

In addition to the preset programs and frequencies (Multi, Alpha, Nogier, Bahr, Reininger), individual patient programs can be easily entered, stored under the patient's name and retrieved. In this manner, you can incidentally prepare your own database. Power (mW), frequencies (Hz), treatment time (min) and number of laser diodes can be varied.

Flexibility when you need it

Flexible treatment arm and intelligent cable routing. Adhesive applicators, which simplify affixing the waveguide and save the



Offer your patients a pain free acupuncture treatment:

- + Simple to use touch screen
- + Simultaneously treat up to 12 acupuncture points
- + Space-saving thanks to compact structure
- + Low life cycle costs



therapist valuable time, are provided as standard. Red laser light (660 nm) is provided as a proven standard and is also available in 785 nm.

Product Specifications:

- Classification Class IIa according to EU Guideline RL 93 / 42 EC
- Laser class 3R according to EN 60825-1:2008
- Effective output in CW operation 12x50 mW
- Number of lasers: 12, up to 12 can be used simultaneously
- Wavelengths: 660 nm or 785 nm available
- Conformity CE 0366
- Pays for itself after only a few treatments
- Made in Germany: guaranteed highest product quality and safety
- Multilingual Product Manual (EN, FR, ES, IT, PT, DE)

• **3B LASER, 12 channel red (660 nm)**
1017713 \$10,488.00

• **3B LASER, 12 channel infrared (785 nm)**
1017714 \$10,488.00



Laser Shower Power Twin 21

Innovation for your practice. The new Power Twin 21 laser shower – holder of the reddot design award – has been specially designed for use on large areas and now offers new treatment options.

For the first time, the laser shower and the laser comb have been brought together in a single device:

- More laser power thanks to 21 laser diodes
- Large beam area of 55 cm²
- Individually adjustable settings via the colour display
- Switch to a laser comb at the click of a button
- Very hairy skin can be treated
- Continuous ray and multi-frequency mode: Nogier, Bahr, Reiningger frequencies have been pre-set
- Laser power can be set at between 5 mW and 50mW
- Available in 1,050 mW (21x50 mW) or 2100 mw (21x100 mw)

• **Laser Shower Power Twin 21; 21 x 50 mw**
1004911 \$4,795.00

• **Laser Shower Power twin 21; 21 x 100 mw**
1018742 \$5,900.00



➤ 3B LASER PENS – Now even more powerful!

Available in 200 mW or 500 mW (both 808nm, infrared). The 500 mW it is the most powerful **3B LASER**. Treatment times are shortened, particularly in frequency mode. In addition to the pre-set frequencies (Nogier, Bahr, Reiningger, Chakra), you can program frequencies quickly from 0.1 to 10,000 Hz, and then store these settings (20 storage locations).

Includes two pairs of laser goggles, two rechargeable battery packs (with a total of 4 batteries), an additional charger and aluminum case.



Product Characteristics:

- Indication of treatment time, laser output, dose, frequency and charge level
 - Memory function for most recently stored values
 - Integrated point-finder
 - 808nm wavelength (infrared laser light)
 - Continuous wave CW
 - Adjustment option for treatment time and dose
- **3B LASER PEN 200 mW, infrared (808 nm)**
1019665 \$3,950.00
- **3B LASER PEN 500 mW, infrared (808 nm)**
1019664 \$4,950.00

MASSAGE TABLES & BOLSTERS



3B Basic Massage Table

All basic models of our massage tables are equipped with a fixed headrest and come with a standard carrying case. Very comfortable PVC vinyl covering, 2 in solid foam padding, oil and water resistant.

Color	Item No.	Price
Blue	W60601B	\$195.00
Black	W60601BK	\$195.00
Burgundy	W60601BG	\$195.00
Green	W60601G	\$195.00



3B Deluxe Stationary Table with Lift Back

The 3B Deluxe Stationary Table with Lift Back is complete with adjustable face rest, storage shelf, and rounded corners.

Color	Item No.	Price
Blue	W60637BL	\$535.00
Black	W60637	\$535.00

3B Basic Stationary Table

The 3B Stationary Table is complete with rounded corners, storage shelf, and arm hammock.

Color	Item No.	Price
Blue	W60636BL	\$495.00
Black	W60636	\$495.00



FEATURES:

- + The Jumbo Full Bolsters are especially good for larger clients, who may require additional support
- + The 3B Full Round Bolster is one of the most popular on the market. Its size makes it versatile and a mainstay for many massage therapists
- + The 3B Comfort Bolster has cutouts making lying facedown more comfortable for many women and making this bolster especially good for pregnancy




3B Jumbo Round Bolster, 8.7 x 25.6 in

Color	Item No.	Price
Blue	W60619JB	\$31.00
Black	W60619JBK	\$31.00
Burgundy	W60619JBG	\$31.00
Green	W60619JG	\$31.00



3B Full Round Bolster, 4.7 x 24.8 in

Color	Item No.	Price
Blue	W60609B	\$25.00
Black	W60609BK	\$25.00
Burgundy	W60609BG	\$25.00
Green	W60609G	\$25.00



3B Comfort Bolsters, 8.7 x 25.6 in

Color	Item No.	Price
Blue	W60623CB	\$50.00
Black	W60623CBK	\$50.00
Burgundy	W60623CBG	\$50.00
Green	W60623CG	\$50.00



All 3B bolsters are covered in a durable vinyl cover that easily wipes clean. Each has a carrying strap at one end for convenient travel.

MASSAGE OILS



➤ Go to 3bscientific.com for more massage oils and accessories



Narayan Balm

These deeply penetrating pure oils relieve sore muscles, aches, pains and stress.

Size	Item No.	Price
Regular strength, 1.5 oz	W67367NBD-1	\$7.99
Regular strength, 1.5 oz - 6 pk	W67367NBD	\$56.34
Extra strength, 1.5 oz	W67367NBX-1	\$7.99
Extra strength, 1.5 oz - 6 pk	W67367NBX	\$56.34



Soothing Touch Narayan Oil

Narayan (nah-RYE-yun) has been used for centuries to soothe muscles. This Sesame oil based blend of Mint, Eucalyptus, and Clove creates a cooling sensation locally while the aromatherapy energizes and soothes the spirit.

Size	Item No.	Price
1 oz	W67367N1	\$9.99
4 oz	W67367N4	\$33.99
16 oz	W67367N16	\$114.99
1 oz - 6 pack	W67367ND1	\$59.94



A. Bon Vital Original Oil

Grapeseed, Olive Oil, Jojoba, Avocado, Soybean, and Safflower oils make this perfect mid-weight massage oil. The thicker consistency allows for lasting glide. This oil is enriched with Vitamins A, E, C, and Pro Vitamin B5. Water dispersible. Unscented. No nut oils. Paraben Free.

Size	Item No.	Price
8oz	W67406O8	\$9.85
Half Gallon	W67406OH	\$28.95
Gallon	W67406OG	\$36.50
5 Gallon	W67406OF	\$204.95

B. Bon Vital Therapeutic Touch Oil

Paraben Free this oil is enriched with Olive Oil, Jojoba Oil, Grapeseed Oil and Arnica extract This lightweight oil is perfect for many modalities including stone massage and aromatherapy. Water soluble for easy cleanup of linens. Unscented. No nut oils.

Size	Item No.	Price
8oz	W67402O8	\$11.95
Half Gallon	W67402OH	\$40.95
Gallon	W67402OG	\$60.75
5 Gallon	W67402OF	\$272.95

C. Bon Vital Naturale Oil

Blend of jojoba, safflower, apricot and peach oils along with natural preservatives give this Paraben Free oil glide with a soft, silky finish. Not greasy. Water dispersible. Unscented. No nut oils.

Size	Item No.	Price
8oz	W67401O8	\$12.95
Half Gallon	W67401OH	\$41.95
Gallon	W67401OG	\$61.75
5 Gallon	W67401OF	\$278.05

3BTAPE



3BTAPE ELITE

3BTAPE Elite features a new fabric, which is made up of a synergistic blend of cotton and synthetic fiber to combine the best features of both worlds. The material offers uni-directional elasticity making it possible to prolong it but not widen it. The unique **3BTAPE Elite** fabric composition allows better moisture release and provides higher endurance for high performance athletes and therapists using it for sport specific treatments.

3BTAPE Elite rolls are 16 ft x 2 in. in length and are available in Beige, Blue, Black, and Pink.

Color	Item No.	Price
Blue	S-3BTEBE	\$11.95
Beige	S-3BTEBK	\$11.95
Black	S-3BTEBL	\$11.95
Pink	S-3BTEPI	\$11.95



3BTAPE Bulk Rolls, 102 ft x 2 in
Whatever your sport, with **3BTAPES**, you can take care of your health and your performance.

**NOW
50% off**

Color	Item No.	Price
Blue	S-3BTBLNL	\$35.00
Beige	S-3BTBENL	\$35.00
Black	S-3BTBKL	\$35.00
Pink	S-3BTPINL	\$35.00



CALL US TODAY!

Discounts available for bulk purchases for medical practices, sports clubs and associations.

3BTAPE

3BTAPE is made from a cotton rayon fiber that is latex free and is hypoallergenic. **3BTAPES** have a skin-friendly, thermoplastic wave-like acrylic adhesive on one side. The wave like adhesive is used to push sweat and moisture away from the tape. This means that **3BTAPES** are water resistance, air-permeable and do not come off in the shower, pool, or sauna. They can be stretched lengthways by approximately 130-140%.

3BTAPE rolls are 16 ft x 2 in. in length and are available in Beige, Blue, Black, Green, Pink, and Yellow.

Color	Item No.	Price
Blue	S-3BTBLN	\$8.95
Beige	S-3BTBEN	\$8.95
Black	S-3BTBK	\$8.95
Pink	S-3BTPIN	\$8.95
Yellow	1012803	\$8.95
Green	1012804	\$8.95



3BTAPE Coated Kinesiology Scissors

The latest innovation in engineering and craftsmanship. These scissors are made from high carbon stainless steel which keeps the blade razor-sharp for years. Perfectly balanced for precision cutting and easily cuts multiple strips of kinesiology tape at the same time. These Kinesiology Taping Scissors are specially coated with black carbon and a fluorine resin which virtually eliminates adhesive residue.

S-3BSCISSOR \$34.95

TOOLS, HOT & COLD THERAPY



A. Acuforce® 2.5

Acuforce® 2.5 uses its own weight (2.5 pounds) to apply pressure, saving wear and tear on the therapist's arms, hands, and fingers.

W51081 \$47.00

B. Acuforce® 7.0

The perfect tool for muscle stripping, cross-fiber friction, and trigger point therapy. The seven-pound Acuforce® 7.0 saves strain on the therapists arms, hands, and fingers.

W51082 \$85.00

C. Acuforce® 3.0i

The Acuforce® 3.0i is a three-pound tool that applies pressure in addition to the benefits of a large ice surface for long-lasting, benefits.

W51083 \$59.00

D. Acuforce® Massage Star

The Wide and Narrow Points to perform trigger Point Therapy and Reflexology.

W51087 \$21.95

Cryoderm Pain Relieving

Fast-acting, deep-penetrating, long-lasting, clinically proven to relieve pain symptoms from soft tissue injuries, muscle strains, sprained ligaments, neurological conditions, arthritic bone and joint conditions, tendonitis, neuritis, bursitis, neuropathy, neuromas, plantar fasciitis and much more. Contains natural ingredients MSM, arnica, boswellia, eucalyptus, ILEX, menthol, peppermint oil.



+ 100% colorless
+ Dye free

+ Paraben free
+ Gluten free

Cold Therapy

Description	Size	Item No.	Price
Gel, Single Bottle	4oz	W54360G4	\$8.00
Roll-On	3oz	W54362R	\$8.00
Gel, Pump Bottle	16 oz	W54360G16	\$27.00
Gel, Pump Bottle	Gallon	W54360GG	\$143.00
Spray Bottle	4oz	W54361S4	\$8.00
Spray Bottle	16oz	W54361S16	\$27.00
Spray Bottle	32oz	W54361S32	\$48.00
Gel, Single Packs	12CT	W54360GP	\$50.00



Heat Therapy

Description	Size	Item No.	Price
Gel, Single Bottle	4oz	W54363T4	\$8.00
Roll-On	3oz	W54363TR	\$8.00
Gel, Pump Bottle	16 oz	W54363T16	\$27.00

Chattanooga™ HotPac®

The original Hydrocollator® moist heat HotPac™ for effective hot therapy. Since their introduction in 1947, thousands of doctors and therapists have insisted on this classic treatment modality.



Description	Size	Item No.	Price
HotPac, Standard	10 x 12 in	W50042	\$10.00
HotPac, Neck	24 in	W50043	\$12.00
HotPac, Oversize	15 x 24 in	W50041	\$22.00
HotPac, Knee/Shoulder	10 x 20 in	W49883	\$19.00

Chattanooga™ ColPac®

There's a ColPac® for every application in cold therapy! Designed especially for heavy clinic use in cold therapy, choose from seven blue vinyl sizes and four black polyurethane sizes. Delivers up to 30 minutes of cold therapy soothing relief from acute pain, swelling and fever. Latex free.



Description	Size	Item No.	Price
Blue Vinyl, Neck	23 in	W50062	\$13.00
Blue Vinyl, Eye	7.5 in	W50065	\$7.00
Blue Vinyl, Quarter Size	5.5 x 7.5 in	W50066	\$6.00
Black Polyurethane, Half Size	6.5 x 11 in	W50070	\$8.00

➤ More sizes available online at 3bscientific.com

Light up your practice with **3B LASERS**

See page 6 to learn more!

The innovative **3B LASER** is classified as low-level laser therapy and is painless and non-invasive. Each laser diode is placed on a specific therapy point, and it radiates with high precision into the tissue which requires treatment.

- Mobile touch screen (control touch pad)
- Side-effect-free treatment
- Broad therapy spectrum
- Made in Germany

Available Options:

3B LASER, Red Laser 660 nm: 1017713

3B LASER, Infrared Laser 785 nm: 1017714



Attractive design and high performance – this is how the new flagship of the **3B LASER PEN** series can be described. Easy to operate, portable, powerful and allows for shortened treatment times, particularly in frequency mode.

- Integrated point finder
- Programmable frequencies
- Includes all accessories for use

Available Options:

3B LASER PEN 200 mW, infrared 808 nm: 1019665

3B LASER PEN 500 mW, infrared 808 nm: 1019664



Order online with code **ACU704***

*Use code ACU704 to get 15% off your online order. Valid one per customer. Offer expires on October 31, 2017.



3B Scientific • 2189 Flintstone Dr. Suite O • Tucker, GA 30084 • 3bscientific.com

Phone: 888-326-6335 • **Fax:** 770-492-0111 • **E-mail:** sales@a3bs.com