







# E-SCOPE<sup>®</sup> Clinical Model



## The **E-Scope (Electronic Stethoscope)**

Clinical Model allows medical professionals to auscultate with increased volume by amplifying sounds up to 30 times louder than an acoustic scope. The E-Scope is popular around the world in clinics, hospitals, research laboratories, and telemedicine settings. Additionally, the E-Scope can be paired with an accessory cord to broadcast sounds through speakers or other audio outputs.

-  One-handed control over volume adjustments
-  Specialist adult diaphragm
-  Automatic shut off after 1.5 to 2 minutes
-  Audio output jack
-  Battery operated (single AAA battery)
-  Supplied with accessory pack

**Cardionics**  
™ The Heart of Auscultation 

Phone: 281-488-5901, 800-364-5901 (USA)  
Fax: 281-488-3195  
Email: [info@cardionics.com](mailto:info@cardionics.com)  
Web: [www.cardionics.com](http://www.cardionics.com)  
Address: 910 Bay Star Blvd, Webster, TX 77598 USA

*E-Scope Catalog #718-7700 (Clinical model)  
#718-7710 (Hearing impaired model)  
Technical Support: Unlimited  
Warranty: One year – parts and labor*