



C-CELIA FEATURE COMPARISON

Enhance emergency obstetric training with the C-Celia simulation system for Cesarean deliveries. Please refer to the following table to choose the correct configuration for your training requirements.



C-CELIA: COMPLETE OBSTETRIC TRAINING SYSTEM FOR CESAREAN DELIVERIES

	Complete Cesarean Delivery Simulator	C-Section Fetal Extraction Trainer	Postpartum Hemorrhage (PPH) Control Trainer	Cesarean Emergency Hysterectomy Trainer
Product number	1022567	1022559	1022561	1022564
C-section delivery stage	Pregnant abdomen, pre C-section	Open abdomen with transverse Cesarean incision	Open, post C-section delivery abdomen	Open, post C-section delivery abdomen
Obstetric features	<p>Emergency Cesarean delivery with:</p> <ul style="list-style-type: none"> • Full transverse (Pfannenstiel incision) C-section • Vertical emergency C-section (running from the umbilicus to the symphysis pubis) <p>C-Section through abdominal wall and uterus with incisions through:</p> <ul style="list-style-type: none"> • Abdominal wall • Muscle • Fascia • Linea alba • Uterus • Amniotic sack 	<p>Difficult Cesarean delivery maneuvers:</p> <ul style="list-style-type: none"> • Manual delivery utilizing fundal pressure • Vacuum assisted extraction • Forceps extraction 	<p>Management of Postpartum Hemorrhage (PPH) after Cesarean delivery:</p> <ul style="list-style-type: none"> • Repair of complex lacerations and uterine artery injuries • Supports PPH control techniques like Bakri Balloon, B-Lynch suture, Hayman suture and O-Leary Stitch • Surgical PPH control: Compressive uterine sutures, uterine artery ligation, hypogastric artery ligation • Non-surgical PPH control: uterine fundal massage, bimanual uterine massage, uterine packing, uterine tamponade, external aortic compression 	<p>Postpartum and Caesarean hysterectomy:</p> <ul style="list-style-type: none"> • Hemorrhage control • Clamping and cutting of uterine art • Fully operable removal of uterus
Functions	<ul style="list-style-type: none"> • Includes twin babies with umbilical cord and placenta • Simulated amniotic fluid 	<ul style="list-style-type: none"> • Includes twin babies with umbilical cord and placenta 	<ul style="list-style-type: none"> • Integrated pressurized bleeding and lubrication system 	<ul style="list-style-type: none"> • Contains uterus, bladder, ureters and detailed representation of uterine arteries
Delivery content	<ul style="list-style-type: none"> • Pregnant abdomen emergency C-section simulator • 2 c-section uteri • 5 abdominal foam inserts • 1 set of twin babies 4.12 lbs (1,87kg) • Repair kit: Silicone glue to allow each uterus to be repaired and reused 5 times (allows 10 total uses) • Amniotic sack (100) • Lubricant gel 	<ul style="list-style-type: none"> • Open abdomen C-section simulator and uterus • 1 set of twin babies 4.12 lbs (1,87kg) • Lubricant gel 	<ul style="list-style-type: none"> • Open (post C-section) abdomen simulator • 2 PPH multi use uteri • Lubricant gel • Artificial blood 1l 	<ul style="list-style-type: none"> • Open (post C-section) abdomen simulator • 10 uteri • Lubricant gel • Artificial blood 1 Gallon (3,8l)



MORE INFORMATION AVAILABLE AT 3BSCIENTIFIC.COM



3B Scientific GmbH • Ludwig-Erhard-Straße 20 • 20459 Hamburg • Germany • 3bscientific.com

Phone: +49 40 73966-0 • Fax: +49 40 73966-100 • E-mail: info@3bscientific.com